



Trigger Nebulizer Cap

Takeing Control of Your Health
with Triggered Bottle

APEX Trigger Nebulizer Cap can deliver small and consistent particle size for targeting delivery to the lungs while wasting less medication and low contamination. With APEX nebulizer, patients of all ages are assured the most efficient aerosol therapy and an improved quality of life.

Care for a Healthy Life

APEX

Trigger Nebulizer Cap

Takeing Control of Your Health with Triggered Bottle

Features ▼



Fine Particle Size(MMAD)

Optimal nebulizer offers fine mist at $<3\mu\text{m}$ (MMAD) with evenly distributed particle sizes.




Ideal for Children

Virtually silent operation and continuous vaporization at various angles.



Low medication waste

Triggered bottle with controller can significantly reduce the environmental aerosol emission to 4%.



Breath enhance technology (one-way valve)



Leak free



Ergonomic bottle design



Lock on – continuing mode

Lock off – discretional operation

Specification

Exhalation resistance	$pep \leq 10\text{-}20\text{cmH}_2\text{O}$
Nebulizing method	Triggered or Continuing mode
Fill Lines	0.5ml, 1ml, 2ml, 3ml, 4ml, 5ml, 6ml
Operating flow rate	7-8 lpm, 50 psi (for Hospital) or with portable compressor (for Homecare)
Air containment	Yes, by the trigger
Latex free	Yes
Flow rate	0.45ml/min
Operating angle	Up to 45 degrees
Fine Particle Fraction (% < 4.8 microns)	84%
Mass Median Aerodynamic Diameter (MMAD)	1.6 μm (CI)/3.9 μm (SLB)
Environmental Aerosol Emission	4%

Distributed by

Apex Medical Corp.
 No. 9, Min Sheng St., Tu-Cheng,
 New Taipei City, 23679, Taiwan
 T +886 2 2268 3100
 F +886 2 2268 6525
 marketing@apexmedicalcorp.com
 www.apexmedicalcorp.com
 Print-2015/04 / All rights reserved
 Specifications may change without notice.